



Product Datasheet

CD261 antibody (orb154377)

Description	Mouse monoclonal antibody to CD261 which is a receptor for the cytotoxic ligand TNFSF10/TRAIL. The a...
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP
Immunogen	Fusion protein containing the extracellular part of TRAIL-R1 and the constant part of the heavy chain of the human IgG1.
Target	CD261
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	O00220
Entrez	8797
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Flow cytometry analysis of TRAIL-R1 expr...



Immunocytochemistry staining (confocal m...