



## Product Datasheet

### CD209 antibody (orb154381)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to CD209 which is a member of the chromosome 19 C-type lectin family that ... |
| <b>Reactivity</b>          | Human                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                            |
| <b>Tested Applications</b> | FC, ICC, IHC                                                                                            |
| <b>Immunogen</b>           | CD209-His-tagged fusion protein                                                                         |
| <b>Target</b>              | CD209                                                                                                   |
| <b>Preservatives</b>       | Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide                                             |
| <b>Concentration</b>       | 1 mg/ml                                                                                                 |
| <b>Storage</b>             | Store at 2-8°C. Do not freeze.                                                                          |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Application notes</b>   | Flow cytometry: Recommended dilution: 1-4 µg/ml.                                                        |
| <b>Isotype</b>             | Mouse IgG2a                                                                                             |
| <b>Clonality</b>           | Monoclonal                                                                                              |
| <b>Purity</b>              | Purified by protein-A affinity chromatography.                                                          |
| <b>Uniprot ID</b>          | <b>Q9NNX6</b>                                                                                           |
| <b>Entrez</b>              | <b>30835</b>                                                                                            |
| <b>Dilution Range</b>      | Flow cytometry: Recommended dilution: 1-4 µg/ml.                                                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                         |



Separation of human myeloid cells (red-f...



Flow cytometry surface staining pattern ...