



Product Datasheet

CD99 antibody (orb154388)


Description	Mouse monoclonal antibody to CD99 which belongs to the CD99 family. It plays role in T-cell adhesion...
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, MC
Immunogen	Human thymocytes
Target	CD99
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG2a kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P14209
Entrez	4267
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Mass cytometry (surface staining) of hum...



Separation of human lymphocytes (red-fil...



Flow cytometry surface staining pattern ...