



Product Datasheet

CD11b antibody (orb154392)

Description	Mouse monoclonal antibody to CD11b
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, IP
Immunogen	Rheumatoid synovial cells and human monocytes.
Target	CD11b
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P11215
Entrez	3684
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human neutrophil granulocytes...



Flow cytometry surface staining pattern ...