



Product Datasheet

BrdU antibody (orb154393)


Description	Mouse monoclonal antibody to BrdU
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IHC-P
Immunogen	Bromodeoxyuridine conjugated to BSA
Target	BrdU
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.
Expiration Date	12 months from date of receipt.



Flow
cytometry
intracellular
staining pa...



Flow
cytometry
multicolor
intracellular
...



Separation of
human BrdU
positive cells
...