

Product Datasheet

TNF alpha antibody (orb154395)

Description

Mouse monoclonal antibody to TNF alpha

Reactivity

Human, Porcine, Primate

Conjugation

Unconjugated

Tested

ELISA, FC, ICC, IHC-Fr

Applications
Immunogen

Recombinant human TNF-alpha

Target

TNF alpha

Preservatives

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Concentration

1 mg/ml

Storage

Store at 2-8°C. Do not freeze.

Note

For research use only

Application notes

ELISA: Biotinylated MAb11 can be used as a detection antibody in combination with capture antibody MAb1. Immunohistochemistry (frozen sections): paraformaldehyde-fixed, saponin-treated frozen tissue sections. Flow cytometry: Intracellular staining.

Isotype

Mouse IgG1 kappa

Clonality

Monoclonal

Purity

Purified by protein-A affinity chromatography.

Uniprot ID

P01375

Entrez

7124

Dilution Range

ELISA: Biotinylated MAb11 can be used as a detection antibody in combination with capture antibody MAb1. Immunohistochemistry (frozen sections): paraformaldehyde-fixed, saponin-treated frozen tissue sections. Flow cytometry: Intracellular staining.

Expiration Date

12 months from date of receipt.