



Product Datasheet

CD148 antibody (orb154402)

Description	Mouse monoclonal antibody to CD148 which is also known as DEP-1, HPTP-eta, and PTP receptor type J (...)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human recombinant CD148 (amino acids 1-444)
Target	CD148
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Isotype	Mouse IgG2b
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q12913
Entrez	5795
Dilution Range	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human neutrophil granulocytes...



Flow cytometry surface staining pattern ...