



Product Datasheet

SCIMP antibody (orb154406)

Description	Mouse monoclonal antibody to SCIMP which also known as Nvl, is a palmitoylated transmembrane adaptor...
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP, WB
Immunogen	Recombinant intracellular part of human SCIMP
Target	SCIMP
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Intracellular staining.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q6UWF3
Entrez	388325
Dilution Range	Flow cytometry: Intracellular staining. Western blotting: Recommended dilution: 1-2 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(intracellular
S...



Western
blotting
analysis of
human
SCIMP...