

Product Datasheet

IFN-gamma antibody (orb154407)

Description	Mouse monoclonal antibody to IFN-gamma which belongs to the type II (or gamma) interferon family. IF...
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP, RIA
Immunogen	Recombinant human interferon gamma
Target	IFN gamma
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	ELISA: Capture antibody in combination with detection antibody 4S.B3.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P01579
Entrez	3458
Dilution Range	ELISA: Capture antibody in combination with detection antibody 4S.B3.
Expiration Date	12 months from date of receipt.