



Product Datasheet

CD195 antibody (orb154409)

Description	Mouse monoclonal antibody to CD195 which belongs to the G-protein coupled receptor 1 family. It also...
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IP, WB
Immunogen	CCR5 peptide (Met1-Lys22) KLH conjugate
Target	CD195
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Western blotting: Do not boil cell lysates prior to SDS-PAGE. Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P51681
Entrez	1234
Dilution Range	Western blotting: Do not boil cell lysates prior to SDS-PAGE. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human CD195 positive lymph...



Flow cytometry surface staining pattern ...