


Product Datasheet

PDE8a antibody (orb154411)

Description	Mouse monoclonal antibody to PDE8a which is a member of the cyclic nucleotide phosphodiesterase fami...
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human PDE8a
Target	PDE8A
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Western blotting: Recommended dilution: 1-2 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	O60658
Entrez	5151
Dilution Range	Western blotting: Recommended dilution: 1-2 µg/ml.
Expiration Date	12 months from date of receipt.



Western blotting analysis of human PDE8A ...