



Product Datasheet

CD85j antibody (orb154415)


Description	Mouse monoclonal antibody to CD85j which is a member of the leukocyte immunoglobulin-like receptor (...)
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	Hairy cell leukaemia cells
Target	CD85j
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG2b kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q8NHL6
Entrez	10859
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human CD85j
positive B
cel...



Flow
cytometry
surface
staining
pattern ...