



Product Datasheet

CD13 antibody (orb154419)

Description	Mouse monoclonal antibody to CD13
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FA, FC, IHC-Fr, IP
Immunogen	Human AML cells
Target	CD13
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P15144
Entrez	290
Dilution Range	Functional application: The antibody WM15 inhibits infection of cells by human coronavirus and inhibits aminopeptidase N activity of the CD13 molecule immunoprecipitates. Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human neutrophil granulocytes...



Flow cytometry surface staining pattern ...