



Product Datasheet

CD11b antibody (orb154428)

Description	Mouse monoclonal antibody to CD11b
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FA, FC, ICC, IHC-Fr, IP
Immunogen	Rheumatoid synovial cells and human monocytes.
Target	CD11b
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P11215
Entrez	3684
Dilution Range	Functional application: Blocking of granulocyte activation. Flow cytometry: Recommended dilution: 1-12 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human neutrophil granulocytes...



Flow cytometry surface staining pattern ...