

Product Datasheet

CD203c antibody (APC) (orb154435)

Description

Mouse monoclonal antibody conjugated to APC which

Reactivity

Human

Conjugation

APC

Tested

FC

Applications

Immunogen

HEK-293 cells transfected with human CD203c

Target

CD203c

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Extracellular and intracellular staining.

Isotype

Mouse IgG1

Clonality

Monoclonal

Purity

Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

[O14638](#)

Entrez


[5169](#)

Dilution Range


Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Extracellular and intracellular staining.

Expiration Date

12 months from date of receipt.



Flow cytometry analysis (surface stainin...



Flow cytometry multicolor staining patte...