

Product Datasheet

BCL2 antibody (APC) (orb154452)

Description

Mouse monoclonal antibody against BCL2 conjugated to

Reactivity

Human

Conjugation

APC

Tested

FC

Applications
Immunogen

Synthetic peptide corresponding to the amino acids 41-54 of human Bcl2

Target

BCL2

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Intracellular staining.

Isotype

Mouse IgG1

Clonality

Monoclonal

Purity

Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

P10415

Entrez


596

Dilution Range


Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Intracellular staining.

Expiration Date

12 months from date of receipt.



Separation of human lymphocytes (red-fil...



Flow cytometry intracellular staining pa...