



Product Datasheet

CD63 antibody (Biotin) (orb154462)


Description	Mouse monoclonal antibody conjugated to Biotin which recognizes is a CD63 that is lysosomal membrane...
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC, ICC, IHC-P, IP
Immunogen	HPB-ALL T cell line
Target	CD63
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-3 µg/ml; positive material: activated platelets, neutrophils and basophils.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P08962
Entrez	967
Dilution Range	Flow cytometry: Recommended dilution: 1-3 µg/ml; positive material: activated platelets, neutrophils and basophils.
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
CD63
positive
CD203c
posit...



Flow
cytometry
surface
staining
pattern ...