


Product Datasheet

CD13 antibody (Biotin) (orb154463)

Description	Mouse monoclonal antibody against CD13 conjugated to Biotin.
Reactivity	Human, Primate
Conjugation	Biotin
Tested Applications	FC, IHC-Fr, IP
Immunogen	Human AML cells
Target	CD13
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Biotinylated antibody is designed for indirect immunofluorescence analysis by Flow Cytometry
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P15144
Entrez	290
Dilution Range	Flow cytometry: Recommended dilution: 3-5 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...