


Product Datasheet

TCR alpha/beta antibody (FITC) (orb154478)

Description	Mouse monoclonal antibody conjugated to FITC which recognizes TCR alpha/beta that recognizes foreign...
Reactivity	Primate, Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	Rat T blasts and erythrocytes
Target	TCR alpha/beta
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-3 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Dilution Range	Flow cytometry: Recommended dilution: 1-3 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...