



Product Datasheet

CD8b antibody (FITC) (orb154480)

Description	Mouse monoclonal antibody conjugated to FITC which recognizes CD8B that also known as P37 and LEU2. ...
Reactivity	Rat
Conjugation	FITC
Tested Applications	FC
Immunogen	CD8 positive Wistar rat splenic T cell hybridomas
Target	CD8b
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P10966
Entrez	926
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of rat CD8b positive cells (r...



Flow cytometry surface staining pattern ...