

## Product Datasheet

### CD99 antibody (FITC) (orb154482)

## Description

Mouse monoclonal antibody conjugated to FITC which

### Reactivity

Human

### Conjugation

FITC

### Tested Applications

FC

### Immunogen

Human thymocytes

### Target

CD99

### Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

### Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

### Note

For research use only

### Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10<sup>6</sup> cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.

### Isotype

Mouse IgG2a kappa

### Clonality

Monoclonal

### Purity

Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

### Uniprot ID

[P14209](#)

### Entrez


[4267](#)

### Dilution Range


Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 10<sup>6</sup> cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.

### Expiration Date

12 months from date of receipt.



Separation of human CD99 positive lympho...



Flow cytometry surface staining pattern ...