


Product Datasheet

CD103 antibody (FITC) (orb154486)

Description

Mouse monoclonal antibody against CD103 conjugated



Reactivity

Human, Primate

Conjugation

FITC

Tested Applications

FC

Immunogen

HTLV-1 induced human T cell line MAPS16

Target

CD103

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.

Isotype

Mouse IgG1 kappa

Clonality

Monoclonal

Purity

Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

[P38570](#)

Entrez

[3682](#)

Dilution Range

Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.

Expiration Date

12 months from date of receipt.

Flow
cytometry
analysis
(surface
stainin...