


Product Datasheet

CD271 antibody (PE) (orb154501)

Description

Mouse monoclonal antibody conjugated to PE which



Separation of
REH cells
stained using
an...

Reactivity

Feline, Ferret, Human, Primate, Rabbit

Conjugation

PE

Tested

FC

Applications

Immunogen

Purified CD271 protein isolated from human melanoma cell line A875

Target

CD271

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Extracellular and intracellular staining.

Isotype

Mouse IgG1

Clonality

Monoclonal

Purity

Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

P08138

Entrez

4804

Dilution Range

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Extracellular and intracellular staining.

Expiration Date

12 months from date of receipt.