



Product Datasheet

CD28 antibody (PE) (orb154508)

Description	Mouse monoclonal antibody conjugated to PE which recognizes CD28 that is expressed constitutively on...
Reactivity	Rat
Conjugation	PE
Tested Applications	FC
Immunogen	Mouse A20J B lymphoma cells transfected with rat CD28
Target	CD28
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P31042
Entrez	25660
Dilution Range	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of
rat CD28
positive
splenocy...



Flow
cytometry
surface
staining
pattern ...