


Product Datasheet

ILF3 antibody (orb154553)

Description	Goat polyclonal antibody to ILF3
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	ILF3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: In paraffin embedded Human Kidney shows nuclear staining in convoluted tubules. Recommended concentration, 2-4 µg/ml.
Protein Sequence	QKPSYGSYGQSHQGG
Clonality	Polyclonal
MW	95.3; 76.0; 74.6
Entrez	3609
Dilution Range	ELISA: 1:32000, IHC-P: 2-4µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Ki...