


Product Datasheet

JAG1 antibody (orb154559)

Description	Goat polyclonal antibody to JAG1, which widely expressed in adult and fetal tissues. It may be invol...
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Target	JAG1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: In paraffin embedded Human Kidney shows membrane staining in the epithelial cells of convoluted tubules. Recommended concentration, 2-4 µg/ml.
Protein Sequence	TNKQDNRDLESAQS
Clonality	Polyclonal
MW	134
Entrez	182
Dilution Range	ELISA: 1:8000, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Ki...