



Product Datasheet

GLP1R antibody (orb154561)


Description	Goat polyclonal antibody to GLP1R, which belongs to the G-protein coupled receptor 2 family. GLP-1R ...
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Target	GLP1R
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 55kDa band observed in Human Brain (Cerebral and Frontal Cortex, Hippocampus) lysates (calculated MW of 53.0kDa according to NP_002053.3). Recommended concentration: 1-3 µg/ml.
Protein Sequence	QHQWDGLLSYQD
Clonality	Polyclonal
MW	53
Entrez	2740
Dilution Range	ELISA: 1:128000, WB: 0.1-1 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Hippocampus lysates.



WB analysis of Rat (A) and Pig (B) Brain...



Immunofluorescence analysis of HeLa cell...