


Product Datasheet

SLC26A1 antibody (orb154562)

Description	Goat polyclonal antibody to SLC26A1 which is a member of sulfate/anion transporter family. This gene...
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	SLC26A1 / SAT-1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: In paraffin embedded Human Kidney shows strong basal membrane staining in convoluted tubules. Recommended concentration, 2-4 µg/ml.
Protein Sequence	QTARARHRELEATD
Clonality	Polyclonal
MW	75.0; 25.0
Entrez	10861
Dilution Range	ELISA: 1:4000, IHC-P: 2-4µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Ki...