


Product Datasheet

FOXF1 antibody (orb154563)

Description	Goat polyclonal antibody to FOXF1, which is a member of the forkhead family. FOXF1 previously design...
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FOXF1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 37kDa band observed in Mouse Lung lysates (calculated MW of 40.0kDa according to Mouse NP_034556.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	HQQVTYQDIKPC
Clonality	Polyclonal
MW	40.1
Entrez	2294
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Lung lysa...