



Product Datasheet

NAT1 antibody (orb154571)

Description	Goat polyclonal antibody to NAT1, which is a member of the arylamine N-acetyltransferase family. It ...
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	NAT1 (aa272-286)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 30kDa band observed in Human Erythrocytes lysates (calculated MW of 33.9kDa according to NP_000653.3). Recommended concentration: 1-3 µg/ml.
Protein Sequence	NISLQRKLVPKHGDR
Clonality	Polyclonal
MW	33.9
Entrez	9
Dilution Range	ELISA: 1:128000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Erythrocyte...



Immunohistochemical staining of Human Ur...