


Product Datasheet

Human Gp130/IL6ST ELISA Kit (orb154573)

Description	Human Gp130/IL6ST ELISA Kit
Reactivity	Human
Range	125 pg/ml - 8,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L24-E619
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	103537 MW
Uniprot ID	P40189
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum and plasma is around 190-446 ng/ml.
Expiration Date	Please enquire.



Human
Gp130/IL6ST
PicoKine
ELISA Kit
sta...