


Product Datasheet

Human CCL13/MCP4 ELISA Kit (orb154575)

Description	Human CCL13/MCP4 ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q24-T98
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	10986 MW
Uniprot ID	Q99616
Sensitivity	10 pg/ml
Dilution Range	For serum, dilution ratio is 1:1 and concentration is around 40 pg/ml. For plasma, dilution ratio is 1:1 and concentration is around 100 pg/ml.
Expiration Date	Please enquire.



Human
CCL13/MCP4
PicoKine
ELISA Kit
stan...