


Product Datasheet

Human MMP-12 ELISA Kit (orb154578)

Description	Human MMP-12 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L17-C470
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	54002 MW
Uniprot ID	P39900
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human MMP-12 PicoKine ELISA Kit standard...