


Product Datasheet

Human CCL14/HCC-1 ELISA Kit (orb154583)

Description	Human CCL14/HCC-1 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: T20-N93
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	10678 MW
Uniprot ID	Q16627
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:500, concentration in serum and plasma is around 60 ng/ml.
Expiration Date	Please enquire.



Human CCL14/HCC-1 PicoKine ELISA Kit sta...