


Product Datasheet

Human CCL16/HCC-4 ELISA Kit (orb154587)

Description	Human CCL16/HCC-4 ELISA Kit
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q24-Q120
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	13600 MW
Uniprot ID	O15467
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
CCL16/HCC-
4 PicoKine
ELISA Kit
sta...