


Product Datasheet

Mouse CCL12/MCP5 ELISA Kit (orb154588)

Description	Mouse CCL12/MCP5 ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: G23-G104
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	11659 MW
Uniprot ID	Q62401
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 80 pg/ml.
Expiration Date	Please enquire.



Mouse
CCL12/MCP5
PicoKine
ELISA Kit
stan...