


Product Datasheet

Human TLR1/TIL ELISA Kit (orb154589)

Description	Human TLR1/TIL ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S22-N578
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	90291 MW
Uniprot ID	Q15399
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 156 pg/ml.
Expiration Date	Please enquire.



Human TLR1
PicoKine
ELISA Kit
standard c...