


Product Datasheet

Rat IL-17F/Interleukin-17F ELISA Kit (orb154592)

Description	Rat IL-17F/Interleukin-17F ELISA Kit
Reactivity	Rat
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A19-A153
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	17712 MW
Uniprot ID	Q5BJ95
Sensitivity	10 pg/ml
Dilution Range	For serum, dilution ratio is 1:1 and concentration is 0-722 pg/ml. For plasma, dilution ratio is 1:1 and concentration is 0-204 pg/ml.
Expiration Date	Please enquire.



Rat IL-17F
PicoKine
ELISA Kit
standard
C...