

Product Datasheet

Human Nesfatin-1/NUCB2/Nucleobindin-2 ELISA Kit (orb154594)

Description Human Nesfatin-1/NUCB2/Nucleobindin-2 ELISA Kit

Reactivity Human

Range 31.2 pg/ml - 2,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: E.coli; Immunogen sequence: V25-L106

Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only


Sample Types cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).

MW 50196 MW

Uniprot ID **P80303**

Sensitivity 10 pg/ml

Expiration Date Please enquire.



Human
Nesfatin-1
PicoKine
ELISA Kit
stan...