

Product Datasheet

PPARG (Phospho-T166) polyclonal antibody [KU] (orb1546576)

Description	PPARG (Phospho-T166) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human PPARG around the phosphorylation site of Threonine 166.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 60 kDa
Uniprot ID	P37231
Dilution Range	WB: 1:1000~1:2000
Expiration Date	12 months from date of receipt.