

Product Datasheet

Histone H3 (Phospho-S28) polyclonal antibody [KU] (orb1546581)

Description	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic phosphopeptide derived from human Histone H3 around the phosphorylation site of Serine 28.
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Concentration	1mg/ml
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
MW	~ 19 kDa
Uniprot ID	P68431
Dilution Range	WB: 1:1000~1:2000 IF: 1:50~1:200
Expiration Date	12 months from date of receipt.