

Product Datasheet

Recombinant Mouse CD40/TNFRSF5 Protein, His Tag (orb1550295)

Description	Recombinant Mouse CD40/TNFRSF5 Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 0.1 EU/ug by LAL method.
Target	CD40
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP2194, Recombinant Mouse CD40/TNFRSF5 Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Leu20-Arg193) of Mouse CD40 (Accession #NP_035741.2) fused with an 6xHis Tag at the C-terminus.
Expiration Date	6 months from date of receipt.