

Product Datasheet

Recombinant Mouse B7-H5/Gi24/VISTA Protein, His Tag (orb1550297)

Description	Recombinant Mouse B7-H5/Gi24/VISTA Protein, His Tag
Species/Host	Human
Reactivity	Mouse
Endotoxins	< 0.1 EU/ug by LAL method.
Target	B7H5
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP2192, Recombinant Mouse B7-H5/Gi24/VISTA Protein is produced by HEK293 Cells expression system. The target protein is expressed with sequence (Phe33 - Ala191) of Mouse VISTA (Accession #NP_083008.1) fused with an 6×His Tag at the C-terminus.
Expiration Date	6 months from date of receipt.