

Product Datasheet

Recombinant Human BAFF trimer Protein, His Tag-Flag (orb1550315)

Description	Recombinant Human BAFF trimer Protein, His Tag-Flag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	BAFF trimer
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP2174, Recombinant mouse FcRn Protein is produced by mammalian expression system. The target protein is expressed with sequence (Thr141-Leu285 Trimer) of human BAFF trimer (Accession #Q9Y275) fused with a 6xHis Tag, Flag tag at the N-terminus.
Expiration Date	6 months from date of receipt.