

Product Datasheet

Recombinant Human FGL1 Protein, Fc His Tag (orb1550328)

Description	Recombinant Human FGL1 Protein, Fc His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	FGL1
Note	For research use only
Purity	> 90% by SDS-PAGE.
Source	KMP2161, Recombinant Human APRIL Protein is produced by mammalian expression system. The target protein is expressed with sequence (Asp64-Asn305) of human FGL1 (Accession #Q08830) fused with an Fc tag at the N-terminus.
Expiration Date	6 months from date of receipt.