

Product Datasheet

Recombinant Human TSLP(R127A, R130A) Protein, Avi His Tag (orb1550353)

Description	Recombinant Human TSLP(R127A, R130A) Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	TSLP
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2136, Recombinant Human TSLP(R127A, R130A) Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Tyr29-Gln159;(R127A, R130A)) of Human TSLP(R127A, R130A) fused with His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.