

Product Datasheet

Biotinylated Recombinant Human TRAIL R4 / TNFRSF10D Protein, Avi His Tag (orb1550371)

Description	Biotinylated Recombinant Human TRAIL R4 / TNFRSF10D Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	TRAIL R4
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2118, Biotinylated Recombinant Human TRAIL R4 / TNFRSF10D Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Ala56-His Tag211) of Human TRAIL R4 / TNFRSF10D fused with His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.