

Product Datasheet

Recombinant Human SOST/Sclerostin Protein, Avi His Tag (orb1550393)

Description	Recombinant Human SOST/Sclerostin Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	SOST
Note	For research use only
Purity	> 95% by Tris-Bis PAGE
Source	KMP2096, Recombinant Human SOST/Sclerostin Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Gln24-Tyr213) of Human SOST/Sclerostin fused with a His Tag and Avi tag at the N-terminal.
Expiration Date	6 months from date of receipt.