

Product Datasheet

Recombinant Human FGL1 Protein, Avi & His Tag (orb1550411)

Description	Recombinant Human FGL1 Protein, Avi & His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	< 0.1 EU/ug by LAL method.
Target	FGL1
Note	For research use only
Purity	> 95% by SDS-PAGE.
Source	KMP2078, Recombinant Mouse CD40 Protein is produced by mammalian expression system. The target protein is expressed with sequence (Asp64-Asn305) of human FGL1 (Accession #Q08830) fused with a 6xHis Tag, Avi tag at the N- terminus.
Expiration Date	6 months from date of receipt.