

Product Datasheet

Recombinant Human Integrin ITGAV & ITGB5 Protein, Avi His Tag (orb1550454)

Description	Recombinant Human Integrin ITGAV & ITGB5 Protein, Avi His Tag
Reactivity	Human
Endotoxins	< 1 EU/ug by LAL method.
Target	Integrin ITGAV & ITGB5
Note	For research use only
Purity	> 95% by Tris-Bis PAGE, > 95% by SEC-HPLC
Source	KMP2035, Recombinant Human Integrin ITGAV & ITGB5 Protein is produced by Expi293 expression system. The target protein is expressed with sequence (Phe31-Val992(P06756) acidic tail&Gly24-Asn719(P18084) basic tail) of Human Integrin ITGAV & ITGB5 fused with His Tag and Avi tag at the C-terminal.
Expiration Date	6 months from date of receipt.